


C5 - Water chloride ion (Cl-) dosage - Water conductivity determination

1. Water chloride ion (Cl-) dosage

1.1. Principle :



Potassium chromate $\text{K}_2\text{Cr}_2\text{O}_4$ in demineralized water is a yellow compound.

It is an end reaction indicator : the first quantities of surplus silver Ag^+ ions, give, with potassium chromate, a red precipitate :



Data :

$$M \text{AgNO}_3 = 108 + 14 + 48 = 170 \text{ g / mol}$$

$$C = 2.818 * 10^{-2} \text{ mol AgNO}_3 / \text{L}$$

$$1 \text{ mL AgNO}_3 \text{ poured} = 2.818 * 10^{-5} \text{ mol AgNO}_3 \text{ react with } 2.218 * 10^{-5} \text{ mol Cl}^- = \mathbf{1 \text{ mg Cl}^-}$$

water
Field of application : 0.5 to 15 mg Cl⁻ in 100 mL water test volume (WTV) :

5 to 150 mg Cl⁻ / L

1.2. Procedure

Test water volume : 100 mL (be careful to the amount of chloride)

Add 1 mL of potassium chromate : solution is yellow

Pour, with burette, silver nitrate until appearance of permanent brown - red precipitate

2. Water conductivity determination

2.1. Principle :

Conductivity is the water property to conduce electric current ; it depends on the kind of dissolved ions, their concentration, and temperature.


Reminder :

$$\text{Resistor (ohms } \Omega \text{)} = \text{Voltage} / \text{Intensity (A)}$$

$$G_{\text{conductance (Siemens S)}} = 1 / R$$

Sensor description : it is constituted by 2 electrodes between which a given voltage is applied.

Resistor measure or electric conductance measure of a water volume delimited by these two parallel platinum electrode :



R is the water volume resistor (Ω), S its surface (cm^2), L its width (cm)

The resistivity is ρ : $\rho (\Omega \cdot \text{m}) = R * S / L$

The **conductivity** γ is : $\gamma (\text{S} \cdot \text{cm}^{-1}) = 1 / \rho = 1/R * L / S$

Other unit : $\mu\text{S} \cdot \text{cm}^{-1}$

- ☞ L / S is called *cell constant* (cm^{-1}) : it must be adjusted as indicated in the notice (generally about 1 cm^{-1}).
- ☞ Conductivity depends on temperature : some conductimeters give directly conductivity values *recalculated for a temperature equal to 20°C (reference)*; these conductimeters have temperature sensor ; other indicate conductivity measured without compensation (test water temperature)
- ☞ Demineralized water conductivity is about $1 \mu\text{S} / \text{cm}$

☞ Standard solutions (20°C) :
 KCl $0.01 \text{ mol} / \text{L}$: $1,277 \mu\text{S} / \text{cm}$
 KCl $0.1 \text{ mol} / \text{L}$: $11,680 \mu\text{S} / \text{cm}$

The recent conductimeters don't need to adjust the zero during calibration.

☞ drinking water conductivity in between about 300 and $600 \mu\text{S} / \text{cm}$; its chloride concentration should be inferior to $250 \text{ mg} / \text{L}$ and its conductivity inferior to $2,500 \mu\text{S} / \text{cm}$.

☞ MM Na = $23 \text{ g} / \text{mol}$
 MM Cl = $35.5 \text{ g} / \text{mol}$

2.2. Procedure

Read the notice conductimeter :

- is there a temperature compensation ?
- calibrate the conductimeter.

Immerse the electrode ; don't stir (conductivity depends on agitation)

Select units and range

Read the conductivity value.

Field of application : drinking waters, industrial waters

3. Practical works

3.1. With powdered NaCl and volumetric flasks, prepare with precision, solutions containing chloride concentrations between about 5 and $150 \text{ mg} / \text{L}$: at least five points

Example of solution preparation :

5 to 150 mg Cl / L = 7.53 to 226.05 mg NaCl / L

preparation of solutions : between 10 and 200 mg NaCl / L

One litre of Stock Solution SS : 2 g NaCl / L

<i>Flask A</i> <i>0.01 g NaCl / L</i>	<i>Flask B</i> <i>0.01 g NaCl / L</i>	<i>Flask C</i> <i>0.1 g NaCl / L</i>	<i>Flask D</i> <i>0.15 g NaCl / L</i>	<i>Flask E</i> <i>0.2 g NaCl / L</i> <i>SS / 10</i>
<i>10 mL E in</i> <i>200mL</i>	<i>50 mL E in</i> <i>200mL</i>	<i>100 mL E in 200</i> <i>mL</i>	<i>150 mL E in 200</i> <i>mL</i>	<i>100 mL SS in</i> <i>1000 mL</i>

3.2. Measure each solution chloride concentration

3.3. Measure each solution conductivity

4. Report

4.1. Draw a board presenting, for each prepared solution :

- a) weighed powdered NaCl masses
- b) prepared solutions volumes
- c) poured AgNO₃ volumes
- d) chloride concentration : mg Cl⁻ / L
- e) conductivity : μS / cm

4.2. Interpret part d) with waited values.

4.3. Plot conductivity = f (chloride concentration)

Verification of your work precision : calculate the correlation coefficient : it must be inferior to 99 %

Material and reagents requirements (for 12 students)

6 flasks containing 400 mL AgNO_3 4.791 g/L

Flask for Ag containing liquid waste

6 flasks with 50 mL of $\text{K}_2\text{Cr}_2\text{O}_4$ 50g/L

Powdered NaCl + weighing machine (precise)

Demineralized water

6 conductimeters with their reference solution + notice

6 volumetric flasks of 1L

30 flasks : 100 or 200 mL