


C8 - Water hardness and water alkalinity measure

1. Water hardness

1.1 Principles

Water hardness is the result of calcium and magnesium (Ca^{2+} and Mg^{2+}) presence.
 Hard water is scale-forming (CaCO_3) in distribution network; it lathers with difficulties.
 Soft water is aggressive, corrosive ($\text{Fe}^0 \rightarrow \text{Fe}^{2+}, \text{Fe}^{3+}$), induces rust formation ($\text{Fe}(\text{OH})_3$); it lathers easily.

Calcium carbonate stability :




Temperature , $\text{Ca}^{2+}, \text{Mg}^{2+}, \text{HCO}_3^-$ concentrations and pH	}	increases induce reaction 1 : water become hard
Temperature , $\text{Ca}^{2+}, \text{Mg}^{2+}, \text{HCO}_3^-$ concentrations and pH	}	decreases induce reaction 2 : water become soft

Calcium and magnesium concentration can be expressed with many different units :

- mol_{ca, mg} / L
- equivalent per liter (eq_{ca, mg} / L)
- g_{ca, mg} / l and mg_{ca, mg} / L
- English, Deutsche or French degrees : these units won't be used in the course.*
- mg CaCO_3 / L or ppm CaCO_3

Correspondence between mol / L, eq / L and ppm CaCO_3
 (mol / L) divided by valence is (eq / L)
 (eq / L) multiplied by 50,000 is (ppm CaCO_3)



Example :

$$2.3 \cdot 10^{-3} \text{ mol Ca}^{2+} / \text{L} = 4.6 \cdot 10^{-3} \text{ eq} / \text{L} = 230 \text{ ppm CaCO}_3$$

$$3.1 \cdot 10^{-4} \text{ mol HCO}_3^- / \text{L} = 3.1 \cdot 10^{-4} \text{ eq} / \text{L} = 15.5 \text{ ppm CaCO}_3$$

Furthermore :

$$1^\circ \text{ French} = 0.56^\circ \text{ Deutsche} = 0.7^\circ \text{ English} = 10 \text{ mg} / \text{L CaCO}_3 = 10 \text{ ppm CaCO}_3$$

$$1^\circ \text{ Deutsche} = 1.786^\circ \text{ French} = 1.25^\circ \text{ English} = 17.86 \text{ ppm CaCO}_3$$

$$1^\circ \text{ English} = 1.438^\circ \text{ French} = 0.8^\circ \text{ Deutsche} = 14.38 \text{ ppm CaCO}_3$$

We will express these concentrations with

$\text{mol}_{\text{ca,mg}} / \text{L}$


eq / L

$\text{mg}_{\text{ca,mg}} / \text{L}$

or ppm CaCO_3 .

Then $10 \text{ ppm CaCO}_3 = 10^{-4} \text{ mol}_{\text{ca,mg}} / \text{L} = 0.5 \cdot 10^{-4} \text{ eq}_{\text{ca,mg}} / \text{L}$; for a water containing both Ca and Mg the concentration expressed with mg / L can't be calculated : we don't know the calcium and magnesium proportions .

Hardness scale :




Field of application : drinking water (the recommended allowances are 900 mg Ca and 400 mg Mg per day), industrial water.


1.2. Water hardness measure

Ethylene diamine tetraacetic acid : EDTA is used during the complexometric dosage


EDTA = YH₄ at pH < 10



YH_2^{2-} at pH = 10



and at pH = 10, YH_2^{2-} forms a complex with Ca^{2+} or Mg^{2+} (chelation) :





First, Ca and Mg can react with NET, added in the tested solution (Volume V) :

Free NET is **bleu** : complexed with calcium or magnesium, NET is **violet** : the tested solution is violet

If $n_{\text{ca,mg}}$ is the molar concentration (unknown) in the tested water (volume V)

Then EDTA, which concentration is precisely known (C_{edta}), is added

EDTA has more affinity for Ca and Mg than NET




When the last Ca, Mg quantities linked to NET disappear, a volume V_{edta} has been poured and the solution turns blue (free NET) :

$$n_{\text{ca,mg}} = n_{\text{EDTA}}$$

$$n_{\text{ca,mg}} = C_{\text{ca,mg}} * V_{\text{pe}} = C_{\text{edta}} * V_{\text{edta}}$$

$$C_{\text{ca,mg}} = C_{\text{edta}} * V_{\text{edta}} / V_{\text{pe}}$$



Procedure :

The sample volume must contain a Ca, Mg quantity allowing a V_{edta} inferior to the burette maximum volume (20 or 25 mL).

example : $V = 100$ mL

Add 4 mL pH10 buffer and a bit of NET : solution is violet

Fill the burette with EDTA 0.01mol / L

Pour EDTA until solution turns blue.

Calcul example :

$V_{\text{water sample}} = 100$ mL,

$V_{\text{edta versed}} = 18$ mL

$C_{\text{ca,Mg}} = 18 \cdot 10^{-3} \cdot 0.01 / 0.1 = 18 \cdot 10^{-4}$ mol / L = 180 ppm CaCO_3 : hard water

⇒ pH 10 buffer and NET in distilled water gives a blue solution : it allows to verify quickly the distilled water quality.

⇒ the hardness measured after ebullition and filtration is called permanent hardness ; it consists of Ca and Mg linked with Cl^- , SO_4^{2-} ...all anions but HCO_3^- , whereas scale has precipitated ; the other hardness is called temporary hardness

⇒ the hardness measured at pH 12 (add 10 mL NaOH 2.5mol / L in place of pH 10 buffer) allows to determinate only the Ca concentration (at pH 12, $\text{Mg}(\text{OH})_2$ precipitate).

2. Water alkalinity


2.1. Principle

Many ions (OH^- , CO_3^{2-} , HCO_3^- ...) can react with H^+ in water (buffer role) but at **neutral pH**, the HCO_3^- (hydrogen carbonate ion) are predominant :

$$\text{HCO}_3^- / \text{CO}_3^{2-} : \text{pKa} = 10.2 \quad \text{Ka} = (\text{CO}_3^{2-}) (\text{H}^+) / (\text{HCO}_3^-) = 10^{-10.2}$$

$$\text{at pH}=7, \quad \text{CO}_3^{2-} / \text{HCO}_3^- = 10^{-10.2} / 10^{-7} = 10^{-3.2} = 6.3 \cdot 10^{-4}$$


HCO_3^- react with acidity (H^+) to give CO_2



The origin of water alkalinity is the presence of HCO_3^- , at neutral pH


HCO_3^- concentration can be measured adding acid in the water sample, with coloured indicator :
 - alkalinity, or methyl orange alkalinity (m-alk), is measured as adding HCl (hydrochloric acid) 0.04 mol / L in a 100mL tested water sample, with methyl orange (turns yellow to red, at pH = 4.5)

Another kind of alkalinity exists, concerning basic waters (containing CO_3^{2-} and OH^-) :
 - phenolphthalein alkalinity (p-alk) is measured as adding HCl (hydrochloric acid) 0.04 mol / L in a 100mL tested water sample, with phenolphthalein ,(turns rose to colourless, at pH = 8.3)




$$\text{H}_2\text{CO}_3 / \text{HCO}_3^- : \text{pKa} = 6.4 \quad \text{Ka} = (\text{HCO}_3^-) (\text{H}^+) / (\text{H}_2\text{CO}_3) = 10^{-6.4}$$

pH evolution in a basic solution in which Hcl is added, with PP and next methyl orange



Then, in drinking waters, $pH = 7$ and the main reaction is



at the equivalence, $n_{H^+} = n_{HCO_3^-}$

If C_a is the acid concentration and V_a the acid volume versed at the equivalence (about $pH = 4$ and m.o turning yellow to red)

If C_{hco_3} is the hydrogen carbonate concentration un known in the tested water sample and V the sample volume :

$$C_a * V_a = C_{hco_3} * V$$

$$C_{\text{hco}_3} = \text{CaVa} / V$$

Alkalinity is expressed in mol $\text{HCO}_3^- / \text{L}$ or mg CaCO_3 / L (ppm CaCO_3 , as hardness) :
 $10^{-4} \text{ mol HCO}_3^- / \text{L} = 10^{-4} \text{ eq} / \text{L}$ (valence = 1) = 5 ppm CaCO_3

Procedure :

V = 100 mL

Add a bit of methyl orange : solution is yellow

Pour HCl 0.04 mol / L until solution turns red

Example :

Va = 14,3 mL

$C_{\text{hco}_3} = \text{CaVa} / V = 0.04 * 14.3 / 100 = 5.72 * 10^{-3} \text{ mol HCO}_3^- / \text{L} = 5.72 * 10^{-3} \text{ eq} / \text{L} = 286 \text{ ppm CaCO}_3$

m-alk = 286 ppm CaCO_3

Water table, drinking waters must present a high m-alk : they won't be aggressive and present high buffer proprieties : generally, m-alk is about 50 (aggressive water) to 400 ppm CaCO_3 (water with calcium carbonate alkalinity) ; here, m-alk = 286 ppm and is enough to present calcium carbonate stability.

3. Practical work : bottling water analysis verification

3.1. Hardness measurement

Read the bottling water analysis and sample a volume (10 to 100mL) which allows a coherent EDTA versed volume (5 to 25 mL) ; a dilution can be necessary

Add 4 mL pH 10 buffer and NET

Pour EDTA until solution turns blue

Calculate water hardness in ppm CaCO_3 and verify your result (bottle etiquette)

Carry on that experiment again : $\Delta A / A_{\text{average}}$ must be < 0.03 (precision test)

3.2. Alkalinity measurement

Read the bottling water analysis and sample a volume (10 to 100mL) which allows a coherent HCl versed volume (5 to 25 mL) ; a dilution can be necessary

Add methyl orange (4 -5 drops) : solution is yellow

Introduce pH meter

Pour HCl 0.04mol/ L until solution turns red ; read pH at the equivalence point

Calculate m-alk, in ppm CaCO_3 and mg $\text{HCO}_3^- / \text{L}$ and verify your result (bottle etiquette)

Carry on that experiment again : $\Delta A / A_{\text{average}}$ must be < 0.03 (precision test)

3.3. Report

3.3.1. Water hardness

- EDTA poured volumes (mL)
- Test precision interpretation
- hardness determination (ppm CaCO_3)
- interpretation of your analysis in comparison with bottling water etiquette one
- is it a soft or a hard water ?

3.3.2. Water alkalinity

- HCl poured volumes (mL)
- pH equivalence
- Test precision interpretation
- m-alk determination (ppm CaCO_3 and $\text{mg HCO}_3^- / \text{L}$)
- interpretation of your analysis in comparison with bottling water etiquette one
- is this water buffered enough ?

4. Material and reagents requirements (12 students)

4L of bottling water (with indicated chemical composition)

6 flasks with 100 mL of EDTA 0.01mol / L

6 flasks with 100mL of pH10 buffer

dissolve 67.5 g of NH_4Cl in 570 mL of ammonia solution 13 mol / L (25% m/m)

*add 5 g of EDTA $\text{C}_{10}\text{H}_{14}\text{N}_2\text{O}_8\text{Na}_2$ **Mg***

complete to 1000 mL with deionized water

6 flasks with 50 mL of methyl orange 5g/L

6 flasks with 100 mL of HCl 0.04 mol / L

12 pH meter + reference solutions (pH 7 and 4) + notice