



C15- Dosage of chlorine in water : Bunsen method

1. Principle

An excess of iodine ion is added to the sample volume containing chlorine




Solution turns brown



Solution turns colourless

Equivalence :



V is the sample volume (L)

V_T is the volume of sodium thiosulphate poured (mL)

C is the unknown chlorine molar concentration (mol Cl_2 / L)

m is the mass of chlorine in the sample volume (mg Cl_2)

M is the chlorine molecular weight : 71 g / mol

C_T is the sodium thiosulphate concentration (mol / L) = 1 / 35.5 mol / L

$$\text{Then : } 2 * m * 10^{-3} / M = V_T * 10^{-3} * C_T$$

$$2 * m / M = V_T * C_T$$

$$m = V_T * (1/35.5) * 71 / 2$$

$$m_{\text{mg Cl}_2} = V_T \text{ ml thiosulphate}$$

At the equivalence, the weight of chlorine (mg Cl_2) in the sample volume is equal to the thiosulphate volume versed

This method is adapted for high chlorine concentration (> 100 mg Cl_2 / L) such as hypochlorite solution (Bleach...)

2. Procedure

Chose a sample volume which contains between 5 and 20 mg Cl_2 ; this volume generally is 10 mL.

All dilution must be carried out with drinking water, not with deionized water.

Add 20 mL of KI 10% (100g/L) and 10 mL of acetic acid : solution turns brown

Pour, with the burette, sodium thiosulphate 1/35.5 mol / L until solution turns colourless (just before solution turns colourless, add starch : solution turns black and the visual transition is easier)

3. Practical work

Determine the chlorine concentration in the Javel water (Bleach) .

Data : 1°Cl (chlorometric degree) = 3.16 g Cl₂ / L

4. Report

Determination of chlorine concentration in bleach : [Cl₂] =

Explain how you have chosen the sample volume, and / or the dilution of the bleach that you have carried out.

Indicate V and V_T .

Calculate the chlorine concentration in bleach, in g Cl₂/L and in chlorometric degree.

5. Required material and reagents

For one experiment :

250 mL of Bleach (Javel water, 12 or 48 °Cl)

100 mL of iodine potassium 10% (100 g KI / L)


50 mL of pure acetic acid

10 mL of starch

2. Dosage of chlorine in water : Titrimetric method with DPD


1. Principle

DPD : N,N-Diethylphénylene-1,4 diamine is oxidised by chlorine (or free chlorine : Cl₂) :



The sample volume turns red

then DPDox is dosed by Mohr's salt Fe²⁺ :



At the equivalence, the solution turns colourless and :

$$2n_{\text{Cl}_2} = n_{\text{mohr}}$$

$$n_{\text{Cl}_2} = n_{\text{mohr}} / 2$$

If V₁ is the volume of Mohr's salt poured, in mL, at the equivalence

C₁ the concentration of Mohr' salt = 2.8 mmol / L

C_{chlorine} the unknown chlorine concentration in the sample volume V = 100 mL

Then

$$C_{\text{chlorine}} V = C_1 V_1 / 2$$

and $C_{\text{chlorine}} \text{ (mol / L)} = C_1 V_1 / 2V$

$$C_{\text{chlorine}} \text{ (mg/L)} = 2.8 * V_1 * (35,5 * 2) / 2 * 100$$

and $C_{\text{chlorine}} \text{ (mg/L)} = V_1 \text{ (mL)}$

In these conditions, the volume of Mohr's salt poured at the equivalence, in mL, is equal to the chlorine concentration in the sample volume, in mg Cl₂ / L.

This method is adapted for low chlorine concentration (< 25 mg Cl₂ / L, i.e. V₁ < 25 mL) such as drinking water (0.1 to 0.5 mg/L...).

☞ chlorine can be combined, for example with ammonia and other nitrogenous compounds, to give chloramines (NHCl₂ , NH₂Cl , NCl₃) ; in order to dose total chlorine (free chlorine and chloramines) , add KI to the sample :

- chloramines oxidise I⁻ to give I₂ ;
 - I₂ can react with DPD to give a red compound (same reaction between Cl₂ and DPD).
- Then Fe²⁺ react with DPD issued from initial Cl₂ and from I₂.

☞ in order to dose chlorine, there is another method : spectrophotometric method

2. Procedure

In a 250 mL erlenmeyer flask, introduce 100 mL of sample.

Add 5 mL of pH 6.5 buffer and 5 mL of DPD : the solution turns red.

Pour Mohr's salt (burette) until solution turns colourless : V₁ mL.

Total chlorine : add 100 mL of sample, 5 mL of buffer, 5 mL of DPD and 1 g of KI and wait for 2 minutes. Verse Mohr's salt : V₁ correspond to the total chlorine concentration.

3. Practical work

Carry out the determination of the chlorine concentration of the sample noted "Cl₂"

4. Report

5. Required material and reagents

For one experiment :

50 mL of pH 6.5 buffer

(KH ₂ PO ₄	46g
Na ₂ HPO ₄	24g
EDTA,2Na (C ₁₀ H ₁₄ N ₂ O ₈ Na ₂ . 2H ₂ O) 8 g/L	100ml
mercuric chloride HgCl ₂	0.02g
deionized water until	1000mL)

5 g of KI

50 mL of DPD

(* In 250 mL of deionized water, mix 2 mL of sulphuric acid (18 mol / L and 1.84 g/mL) and 25 mL of EDTA,2Na 8g/L

* dissolve in this mixture, 1.1 g of anhydrous diethyl - p - phenylenediamine DPD sulphate [NH₂-C₆H₄-N(C₂H₅)₂ - H₂SO₄]

* complete with deionized water until 1000mL

☞ maximum storage time : one month

50 mL of Mohr's salt $[(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2, 6 \text{H}_2\text{O}]$ $2.8 * 10^{-3} \text{ mol / L}$

$[(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2, 6 \text{H}_2\text{O}]$ 1.106 g

sulphuric acid diluted 1/3 (V/V) 1mL

deionized water until 1000mL)

250 mL of diluted bleach ($0.1 < [\text{Cl}_2] < 25 \text{ mg Cl}_2 / \text{L}$)

For example, if initial chlorine concentration in bleach is 48°Cl (152 g / L) : dilute it 10,000 times (until 15.2 mg / L : 100 μ L of bleach in one litre of drinking water , not in deionized water).


one microburette (0.1 to 1 mL or to 5 mL)

3. Determination of the break point


1. Principle

The break point is the amount of chlorine necessary to remove ammonia in water.


Chlorine first react with ammonia to give chloramines (A) :



Beyond this formation, the chlorine added reacts with chloramines and destroy them (B) ; the amount of chlorine necessary to form and next to remove chloramines is the break point ● ; beyond this point, the amount of added chlorine remains in water as free chlorine (C) .



Evolution of residual chlorine in a sample of water containing 1 mg N.NH₃ / L, in which chlorine is added



7.6 mg of dichlore are necessary to remove 1 mg N. NH₃.

2. Procedure

You have twelve 500 mL erlenmeyer flasks containing 500 mL of solution 1 mg N.NH₄⁺ / L (3.82 g NH₄Cl / L).

Add these amounts of bleach :

flask n°	1	2	3	4	5	6	7	8	9	10	11	12
Bleach 1g/L (mL)	1.5	2	2.5	3	3.5	3.75	4	4.25	4.5	5	5.5	6
final [Cl ₂] mg/L	3	4	5	6	7	7.5	8	8.5	9	10	11	12

Keep the flask in the darkness during 2 hours.

Add 10 mL of acetic acid and 20 mL of KI 10% : solution must turns brown
Titrate total chlorine with sodium thiosulphate 1/35.5 mol / L until solution turns colourless :
note V_T .

3. Report

Plot the curve representing $V_T = f ([Cl]_2)$

Determinate the break point.

4. Material and reagent requirements

12 erlenmeyer flasks : 500 mL , containing 500 mL of 3.82 g NH_4Cl / L

200 mL of acetic acid

300 mL of KI 10%

Starch

1L of sodium thiosulphate 1/35.5 mol / L