

Mi1 : Use of the microscope and wet-mount technique

Requirements for one experiment

One microscope
5 glass slides
5 coverslips
5 mL of broth : E.coli suspension
5 mL of activated sludge suspension
5 non sterile Pasteur pipette
One pipette filler for pipette pasteur

1. Description of the microscope

Page 1 : Comparison of microscope types ; size of cells and microscopic applications.

Page 2 : Bright-field microscope components

2. Wet-mount technique : page 3 (figures 1 to 3)

You have two broth :

- a pure suspension of E.coli, which should not be contaminated : all the experiments must be done aseptically : *observe how to handle tubes and pipettes aseptically* ;
- a suspension of activated sludge : no need to work aseptically.

Carry out the observation of

- a suspension of E.coli (1 to 5 μ m) : with a sterile pipette, draw aseptically one drop of this suspension and place it on a glass slide ; place the edge of the cover slip on the drop ; observe with the objective X40.
- a suspension of activated sludge (5 to 100 μ m) : objective X 10 and then X 40.

After observation, put glass slides and pipettes in bleach and suspensions in the basket destined to sterilisation.

3. Report

For each observation, note :

- the magnification used
- for cell :
 - their average size (μ m)
 - their shape (rods, cocci for bacteria)
 - their association type (single, grapelike clusters, chain...)
 - their mobility

Draw a sketch illustrating your observations.