

**Mi4 : Fecal contamination indicators : thermotolerant and total coliform  
and streptococci in drinking water**

**Filtration technique**

**1. Generalities**

The intestinal flora includes several germs, indicated following this decreasing predominance :

- rod shaped Gram negative bacteria, anaerobic and called bacteroids
- the genus Bifidobacterium ( rod shaped, Gram positive, anaerobe)
- genus Clostridium (cocci, gram positive and sporulating) and lactic bacteria

These three kind of bacteria are difficult to grow in laboratory.

- the specie E.coli
- the genus Enterococcus

These two types of bacteria are easy to grow

And finally,

- enterobacteria (like Citrobacter), genus Staphylococcus, genus Bacillus, yeasts...

FCI are commensally of the intestine ; when they are present in water, there is a risk of contamination due to pathogenic bacteria like Salmonella, Shigella, Vibrio cholerae...

Fecal contamination is due to the discharge of waste water in surface water ; MWW contain  $10^6$  to  $10^8$  bacteria / mL.

There are three categories of FCI :

- total coliforms, thermotolerant coliforms and E.coli :they are commensally of the intestines and can't survive for a long time in water ; E.coli is the main representant of thermotolerant coliforms
- fecal streptococci (Enterococcus faecali) , also commensally of the intestine, but can survive longer in water : they are the indicator of an older contamination
- sulphite reducing Clostridium, which are the less reliable FCI because they can also live like saprophytes in water.

The most reliable FCI in water are thermotolerant coliforms.

**2. Principle and definitions**

**Coliforms :**

Micro-organism growing at 37°C on lactose-bile salt containing agar (called tergitol 7 TTC agar), acid producing in 24h and oxidase negative .

. The selective agar medium Tergitol 7 TTC contains molecules that inhibit cocci gram positive growth (like streptococci...)

**Thermotolerant coliform :**

Same proprieties than coliforms, but at 44°C.

*There are two steps : presumption and then confirmation*

*- growing on tergitol 7 TTC*

*- inoculating each colony from Tergitol 7 on nutritive agar in order to examine the presence of oxidase :*

*We will just carry out the first step, i.e. presumption*

**Streptococci :**

Micro-organism growing at 37°C on glucose- azide containing agar (called Slanetz Bartley agar), generating typical TTC reducing colonies, and generating positive reaction in 24h at 37°C on esculine –bile containing agar. The selective agar medium Slanetz Bartley contains molecules that inhibit rod shaped gram negative growth (like coliforms...)

*We will just carry out the first step*

### 3. Operation

Requirement for one experimentation :

*One flask containing 300 mL of 100 to 300 E.coli / L suspension*

*One flask containing 300 mL of 100 to 300 Enterococcus faecalis / L suspension*

**How to prepare a solution containing a precise bacteria concentration ?**

- inoculate 10mL of broth with the tested bacterium

- incubate 24 h at 37°C :

**the concentration is maximum and  $\approx 10^7$  /ml =  $10^{10}$  /L**

**Dilute this suspension  $10^8$  times : 100  $\mu$ L in 100 mL, then 100 $\mu$ L in 100mL, then 1 in 100mL and you'll have a concentration equal to 100 bacteria /L**

*3 sterile filtration membranes 55mm, 45 $\mu$ m*

*One filtration ramp*

*Two tergitol 7-TTC agar plate 55mm*

*One slanetz barley agar plate 55 mm*

*Two 100 mL flask of sterile distilled water*

*Two incubators : 37 and 44°C*

- Prepare aseptically the filtration apparatus and the membrane : **page 1**

and **page 2** : figures 1 to 5 (but no absorbent pad)

- filtrate 100ml of water sample (3 samples)

- place the membrane : (figures 6 to 9 **page 2** but between figure 7 and 8, rinse the funnel with sterile water)

- on Tergitol 7 TTC agar plate and incubate at 37°C in order to analyse total coliform

- on Tergitol 7 TTC agar plate and incubate at 44°C in order to analyse thermotolerant coliform

- on Slanetz Bartley agar plate and incubate at 37°C in order to analyse fecal streptococci

Duration of incubation : 24h

Second day : Count positive colonies :

- on Tergitol 7 : yellow colonies surrounded by a yellow halo

- on Slanetz : red pink or brown colonies

### 4. Report

Each colony is issued from one bacterium growth

We'll consider that each characteristic colony (first step, i.e. presumption step) represents the definitive analysed bacterium (no confirmation test).

Determine the : (bacteria / 100mL)

- total coliform concentration

- thermotolerant concentration

- fecal streptococci concentration, in the three samples.

*The result should be equal to **approximately** 10 colonies / plate, i.e. 10 bacteria / 100mL. (between 1 and 100 bacteria / 100 mL)*